

Electro Surgical Units (cautery)


ORIS Series

ORIS 400 A



ORIS 400 A is 400W analog model with number of cut, blend and coagulation mode with to do under water surgery. This is the model with all around capabilities and simple operation.

ORIS 400



ORIS 400 is micro controller based model with all the required features. Full power output and low cost. This is right model for the customers who required basic 400w micro controller based model.

ORIS 400 Plus



ORIS 400 Plus is a top of the line micro controller based model with 20 programs, high speed tissue response circuits and host of modes and superior bipolar performance for the customers who demand the best.



Technical Specifications

Features Models	Body	Max RF Power Watt	Model Type	No. of Programmes	Tissue Response Circuit	Patient Plate Monitor	Cut Modes	Coag Modes	Bipolar Mode	Weight	Dimension	Standard Accessories
ORIS 400 Plus	Plastic	400	Micro Controller Based	20	Present	Present	PC+b1+b2+EC	Sp+Dc+Ful	Mic+Mac+Cut	3.8 kg	4.5"x10"x12"	sspp+ch+fs+hpd+bp
ORIS 400	Plastic	400	Micro Controller Based	-NA	Present	Present	PC+b1+b2+EC	Sp+Dc+Ful	Mic+Mac+Cut	3.8 kg	4.5"x10"x12"	sspp+ch+fs+hpd
ORIS 400 A	Plastic	400	Analog	-NA	NA	Present	PC+b1+b2	Sp+Dc+Ful	Single mode	3.8 kg	4.5"x10"x12"	sspp+ch+fs+hpd